

Here is the link to download the free software (NI ELVISmx) needed for the NI myDAQ:

1. Download NI-ELVISmx 19.0
 - a. Browse to <https://www.ni.com/en-us/support/downloads/drivers/download.ni-elvismx.html>
 - b. Select DOWNLOAD
2. Run ni-elvismx_19.0_online_repack.exe
 - a. Accept default "Additional items you may wish to install:" PLUS select to install the Firmware option.
 - b. You will need to restart your computer when it is finished
3. Plug in the NI myDAQ into your computer's USB port.
4. Run Start | NI ELVISmx Instrument Launcher
Select Oscilloscope, Digital Multimeter, Function Generator, etc.

